

LIFEPAK® CR2 AED manual maintenance inspection update LIFELINKcentral™ AED Program Manager

Contents

- Instructions to log an inspection for one AED at a time.
- Instructions to log an inspection for all AEDs at once.

To log an inspection for one AED at a time:

1. Visually confirm your AED(s) flashing green indicator.
2. Log in
 - a. Log into LIFELINKcentral AED Program Manager.
3. Select Site
 - a. From “Sites” select the site for which you would like to log the inspection (Redmond – North Building).

The screenshot displays the LIFELINKcentral AED Program Manager interface. The top navigation bar includes the Physio Control logo, the text "LIFELINKcentral™ AED PROGRAM MANAGER", and a user welcome message "WELCOME Riggs, Malorie". The left sidebar shows a navigation menu with "Alerts" (My Account, Sites), "Equipment", "People", and "Resource Center". The main content area is titled "Sites" and features a search bar, a list of sites, and a "NEEDS ATTENTION (DEVICE)" button. The sites listed are "Redmond-North Building (SS774426)" and "Redmond-South Building (SS776691)".

Site Name	Address	Primary Site Coordinator	Phone Number	E-mail
Redmond-North Building (SS774426)	11811 Willows Road NE Redmond, Washington, 98052	Riggs, Malorie	425-867-4700	mriggs@sample.com
Redmond-South Building (SS776691)	11811 Willows Road NE Redmond, Washington, 98052	Graubard, Angela	425-867-4700	agraubard@sample.com

4. Click the Inspect button to record an inspection.

The screenshot displays the LIFELINKcentral AED Program Manager interface. The top navigation bar includes the Physio Control logo, the title 'LIFELINKcentral™ AED PROGRAM MANAGER', the user account 'Demonstration Customer Account / Redmond-North Building', and a welcome message for 'Riggs, Malorie'. A left sidebar contains a menu with 'Alerts' and 'My Account'. The main content area shows the 'Site' details for 'Redmond-North Building (SS774426)'. An 'Inspect' button is highlighted with a yellow box. Below the site name, a list of site details is provided, including account information, coordinator, phone number, creation date, address, local EMS agency, and AED post-event report recipient.

Site	
Redmond-North Building (SS774426)	
Account	Demonstration Customer Account (AETT0901)
Primary Site Coordinator	Riggs, Malorie
Site Phone Number	425-867-4700
Create Date	7/15/2016
Address	11811 Willows Road NE Redmond Washington 98052 United States
Local EMS Agency	Not Specified
AED Post-Event Report Recipient	Dames, M.D., Joanne
Required number of Trained Responders per Site	1

5. Select AED

- a. Select the AED for which you are recording the inspection.

LIFEPAK CR2 (48248214)

Performed By	System
Readiness Time	Not Specified
Location	First Floor Lobby

NO DATA AVAILABLE

6. Verify information (important for first inspection, optional for future inspections).

- a. Expand each section and verify the information.

Perform Inspection ×

LIFEPAK CR2 (48248214)

EQUIPMENT INFORMATION ▼

ACCESSORIES ▼

INSPECTION CHECKLIST ▲

Device Present ☒

Visible Damage ☐

Readiness Indicator ☒ Flashing ☐ Off

Override Status ☐

Readiness Status Ready

NOTES ▼

Submit **Cancel**

7. Confirm AED readout
 - a. Proceed to Inspection Checklist and confirm the flashing readiness indicator is displayed.

Perform Inspection ✕

LIFEPAK CR2 (48248214)

EQUIPMENT INFORMATION ▼

ACCESSORIES ▼

INSPECTION CHECKLIST ▲

Device Present ☒

Visible Damage ☐

Readiness Indicator ☒ Flashing ☐ Off

Override Status ☐

Readiness Status Ready

NOTES ▼

Submit Cancel

- b. Click the “Submit” button to submit your entry.
 - c. The AED status will turn green for “Ready”, confirming a successfully logged inspection:

LIFEPAK CR2 (48248214)

Performed By Riggs, Malorie

Readiness Time 4/29/2020 2:41 PM

Location First Floor Lobby

READY

You will return to the list view of the AED Inspection Schedule page. From there, repeat the steps above to record a maintenance inspection for the next AED in the list.

To log an inspection for all AEDs at once:

8. Visually confirm the AEDs display the flashing green readiness indicator.
9. Log in
 - a. Log into LIFELINKcentral AED Program Manager.
10. Inspections
 - a. From “Equipment,” click “Manage Equipment”, and then click “Inspections” to see the inspection schedules for those AEDs for which you are responsible.

The screenshot displays the LIFELINKcentral AED Program Manager interface. The top navigation bar includes the Physio Control logo, the title "LIFELINKcentral™ AED PROGRAM MANAGER", and a user welcome message "WELCOME Riggs, Malorie". The left sidebar contains a navigation menu with "Alerts", "My Account", "Sites" (2), "Equipment" (4), "Manage Equipment", "AEDs Ready 3 of 4", "People" (5), and "Resource Center". The "Equipment" section is active, showing a search bar and a dropdown menu with options: "Inspections", "Device Issues", "Custom Equipment", "Accessories", and "AED Usage Log". The main content area displays equipment details for "Redmond-North Building (SS774426)" and "Redmond-South Building (SS776691)".

Equipment

Enter text to search the list below

Redmond-North Building (SS774426)

Equipment ID	Location	Device	Expiration Date
CR2 (48248214)	Lobby	CR2 Lithium Battery	4/29/2022
		CR2 QUIK-STEP Electrodes, Adult/Infant/Child	4/29/2022
CR2 (48248216)	Second Floor Break Room	CR2 Lithium Battery	4/29/2022
		CR2 QUIK-STEP Electrodes, Adult/Infant/Child	4/29/2022

Redmond-South Building (SS776691)

Equipment ID	Location	Device	Expiration Date
LIFEPAK CR2 (47987626)	Warehouse	CR2 Lithium Battery	4/10/2023
		CR2 QUIK-STEP Electrodes, Adult/Infant/Child	10/15/2023

11. Inspection Schedules

- a. Next, click on the inspection schedule to which you would like to navigate. In this example, we will click on “Redmond – South Building,” to log an inspection for all AEDs affiliated with that schedule:

**PHYSIO
CONTROL**

LIFELINKcentral™
AED PROGRAM MANAGER

Demonstration Customer Account / All Sites

WELCOME
Riggs, Malorie

Alerts

My Account

Sites2

Equipment4

Manage Equipment>

Ready	3
No Data Available	1

←

Inspection Schedules

+

Enter text to search the list below

Redmond-North Building (SS774426)

AED Inspection Schedule

Person Responsible	Riggs, Malorie
Frequency	Monthly
Last Inspection Date	4/29/2020

Redmond-South Building (SS776691)

AED Inspection Schedule

Person Responsible	Riggs, Malorie
Frequency	Monthly
Last Inspection Date	4/3/2020

12. AED Inspection Schedule

- a. Once you are on the landing page for the “Redmond – South Building” inspection schedule, click “Inspect All”.

The screenshot displays the LIFELINKcentral AED Program Manager interface. The top navigation bar includes the Physio Control logo, the LIFELINKcentral AED PROGRAM MANAGER title, the user account information (Demonstration Customer Account / Redmond-South Building), and a welcome message (WELCOME Riggs, Malorie).

The left sidebar contains the following menu items:

- Alerts
- My Account
- Sites (2)
- Equipment (2)
- Manage Equipment (>)
- AEDs Ready 2 of 2 (indicated by a green circle)
- Ready (2)
- People (1)

The main content area is titled "AED Inspection Schedule" and features a search bar with the placeholder text "Enter text to search the list below". Below the search bar, there are two AED inspection entries, each with a "READY" status button:

- LIFEPAK CR2 (48125068)**
 - Performed By: System
 - Readiness Time: 4/3/2020 1:57 AM
 - Location: First Floor Lobby
 - Status: READY
- LIFEPAK CR2 (47987626)**
 - Performed By: System
 - Readiness Time: 4/9/2020 8:47 AM
 - Location: Warehouse
 - Status: READY

The "Inspect All" button is highlighted with a yellow border.

13. Confirm

- a. The following prompt will appear, asking you to confirm if you would like to update the readiness status of all AEDs associated with this schedule to “Ready”. Click “Yes” to complete the check:

Complete Inspection

×

The readiness status of all Equipment not previously recorded today will be set to 'Ready' with today's date.
Do you want to proceed?

Yes

No

14. Completed

- a. You will be routed to the list view for all AEDs on that inspection schedule. The maintenance inspection for all AEDs under this schedule is complete, and all AEDs will show green for “Ready” with today’s date.

PHYSIO
CONTROL

LIFELINKcentral™
AED PROGRAM MANAGER

Demonstration Customer Account / Redmond-
South Building

WELCOME
Riggs, Malorie

Alerts

My Account

Sites2

Equipment2

Manage Equipment>

AEDs Ready
2 of 2

Ready2

People1

←AED Inspection ScheduleInspect All↶↷

Enter text to search the list below

LIFEPAK CR2 (48125068)

Performed BySystem

Readiness Time4/3/2020 1:57 AM

LocationFirst Floor Lobby

READY

LIFEPAK CR2 (47987626)

Performed BySystem

Readiness Time4/9/2020 8:47 AM

LocationWarehouse

READY

For further information, please contact your Account Support Representative or call 888-627-0007